

## Fujifilm Drypix Plus

### Description

The DryPix Plus is a dry laser imager which combines the high quality and performance of Fujifilm's centralized imagers, with an extremely compact design and convenient 110 volt power, making it ideal for small centralized departments that do not need the ultra-high speed and throughput of DryPix 5000. The DryPix Plus has excellent throughput between 110 and 160 films per hour. Like the DryPix 5000, this unit also features 50-micron high resolution printing and has been approved for use with Full Field Digital Mammography (FFDM) devices. It accommodates all common film sizes, has three film drawers and an optional top-mounted 4-bin film sorter. An internal DICOM print server with Windows XP operating system allows the DryPix Plus to operate with high performance in environments where space is at a premium.

#### Date

2024/10/06

#### Date Created

2024/05/15

#### Author

sws-dev