



GE Vivid S5

Description

The **GE Vivid S5** is an Ultrasound system for Cardiovascular imaging. The Vivid S5 can be shared with other applications; Musculoskeletal, Abdominal, Pediatric, and Obstetrics. The Vivid S5 provides is known for its speed in producing quality images. The Ultrasound system offers an easy-to-manuever ergonomic design, making for a more efficient, less tree situation. The GE Vivid S5 offers a flexible system. This allows the control panel and the keyboard to be moved around and adjusted when needed. The system offers imaging modes that include 2D Mode, Color Doppler Mode, M-Mode, PW/CW Doppler, and LVO contract. The system offers additional packages for Stress Echo for measuring the stress on a patient's heart.

- Color Doppler Ultrasound
- 2-D, M-Mode, PW/CW Doppler
- Coded Octave Imaging (Harmonics)
- ECG Module & 3-lead cable
- Anatomical M-Mode Color Angio
- EchoPac System post-processing and data archiving
- High-Resolution 15" monitor
- Fully Digital
- Integrated Digital Stress Echo (optional)

Date

2024/10/06

Date Created

2024/05/15

Author

sws-dev