



Konica Sonoimage HS1

Description

The Konica Sonoimage HS1 is a high-end ultrasound system primarily used in medical imaging.

Scan method: Convex, Linear and Sector

Monitor: 15-inch IPS monitor

Dimension: W 369mm x D 452mm x H 90mm (when folding the monitor)

Power consumption: AC100V, 50/60Hz, max. 180VA (when using AC adopter)

Weight: 7.8kg \pm 10% (System only and including battery)

Date

2024/10/06

Date Created

2024/05/15

Author

sws-dev