



## Toshiba Xario 100

### Description

#### 3D/ 4D Toshiba Xario 100 Ultrasound Machine, Color Doppler Velocity

- Model Name/Number. Xario 100.
- Brand. Toshiba.
- Doppler Mode. Color Doppler Velocity.
- Imaging Mode. 3D/ 4D.
- Clinical Application. General Imaging.
- Monitor Size. 19" LCD.
- Type of Transducer/ Probe. Linear Array, Convex, Phased Array.
- Monitor Resolution. 1280X1025.

#### Date

2024/10/06

#### Date Created

2024/05/15

#### Author

sws-dev